



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com



2016 CATALOG

PRESS RELEASE – 1st In the Industry

**Cytotoxicity Tested and Titanium Oxide Free Speculums
and DEHP Free Tubing California Proposition 65 Compliant**



CLEARWATER IS GOING

GREEN




COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com



LABELING

ADEQUATE DIRECTIONS FOR USE:

Under no circumstances should the patient be left unattended.

21CFR801.109

INDICATION FOR USE STATEMENT:

The indication for use of this device must be restricted to colon cleansing when medically indicated, such as before radiological or endoscopic examination.

21CFR876.5220

PRESCRIPTION STATEMENT:

“Caution: Federal law restricts this device to sale by or on the order of a physician.”

21CFR801.109



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON PREP EQUIPMENT AND SUPPLIES

TABLE OF CONTENTS

Clearwater Colon Hydrotherapy Equipment

Labeling	1
Welcome Letter	5
Background	6
PPC-101 Clinical Model with and w/o Specimen Collector	7
PPC-101 Clinical Model with Built In UV Filtration System	8
PPC-101 Clinical Features	9
PPC-101 “The Traveler” Features	13
PPC-101 “The Traveler” Dual Purpose	14
PPC-101 “The Traveler” with Sewage Tank Stand	16
“The Traveler” Stainless Steel Stand	16

Water Filtration System

Granular Activated Carbon Cartridges (10 Inch)	17
(1 Micron Taste/Odor/Chemical) (5 Inch)	17
Fiber Cartridges (10 Inch)	17
(Rust/Sediment 1 Micron (5 Inch)	17
Carbon Filter (Gac Inline) for Other Manufacturer Units (10 Inch) ..	17

Disposable Speculums

Speculum Inflow water Line Disconnects	18
Bio-Compatibility Test Results	19
Speculum Features	19
AutoShip Program	20



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Therapy Table & Accessories

Adjustable Electric Therapy Table	22
Light Weight Portable Tables	23
Table Cart Carrier	23
Pillow Wedge with Vinyl Cover	23
Bolster Round Supports (Knees/Ankles)	23
Therapist's Chair	24
Step Stool	24
Table Cover	24
Table Cover Warmer	24

Colon Therapy Charts

Colon Model	25
Colon Therapy	25
Digestive System	25
Blueprint for Health Chart	25
Diseases of the Digestive System	26
Gastroesophageal Disorders and Digestive Anatomy	26
Bowel Reflex Points	26
Inflammatory Bowel Disease (IBD)	26
Understanding Irritable Bowel Syndrome (IBS)	27



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Products & Supplies

Stainless Steel View Tube Brush Holder 27

Gowns and Footies 27

Books (Complete List in Product Info & Price Section pg 33) 28

Supply Packet 28

Speculum Sample Pack 28

Gloves 28

*Federal Code of Regulations on Purchasing Device and
Speculums*

Sample Letter/Script on Doctors Letterhead 29

Boosting Your Business - Using Video E-Mail

Video E-Mail Program 30

Product Information and Price List

Item Price List 31



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

SOME OF THE REASONS WHY OUR EQUIPMENT WILL GREATLY BENEFIT YOU AND YOUR PATIENTS!

Equipment used for cleansing the colon - such as, before radiological or endoscopic examinations. Used for preparation for a colonoscopy.

- ◆ Designed for use in Gastroenterology Clinic, Hospitals & Nursing Homes
- ◆ Cuts down on bad preps
- ◆ Portable can go to bedside of Colonoscopy table
- ◆ Sewage Tank Stand is great when plumbing can not be done

**Abstract on Colon Cleaning VS Prep Solutions when doing Colonoscopy
(Copy of Abstract upon request)**

Do you have bad preps?

- ◆ Our Equipment monitors the Colon pressure throughout the whole session
- ◆ Completely closed system
- ◆ All Stainless Steel Cabinets
- ◆ Pressure & Temperature Controlled with safety features
- ◆ Water in View Tube is full at all times for viewing and inspecting specimens
- ◆ Safe connects and disconnects of Water, Sewer, and Electrical
- ◆ Options for Water Line and Waste Line connections
- ◆ Quick and easy connector device on each and every Speculum kit line
- ◆ Disposable Speculums that comes in different sizes and styles
- ◆ Speculums are designed for easier insertion and greater client comfort
- ◆ Bio-Compatibility test results study - given Superior Medical Rating
- ◆ Free 2 full days of Training with purchase of a unit

Do you have different locations?

- ◆ Our Equipment is the answer. We have the Patent on Portability
- ◆ Patent on the ONLY Speculum (rectal nozzles) that are medical grade plastic and designed anatomically correct.
- ◆ ONLY inflow water line that is PVC FREE (Important for your sensitive patients)



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Mary Ruth Baker, Clearwater Colon Hydrotherapy owner, has been in the field of Colon Hydrotherapy for over 30 years. Her commitment has always been to have the highest standard of quality and client care. These high standards are now available in Clearwater Colon Hydrotherapy systems. These systems have the best mechanical and operational design, which makes them stand alone as the most efficient colon prep for Colonoscopy, colon cleansing equipment that is on the market today.

Mary Ruth is anxious to teach and support you in the operations of the Clearwater Colon Hydrotherapy equipment and also with the knowledge of the business operations towards setting up a complaint office.

An installation manual & DVD is provided when you purchase a unit, that demonstrates the complete installation process. A two day training seminar on the installation, operation and maintenance of the equipment is provided free of charge, if you schedule your training within 30 days of your initial equipment purchase, at our corporate headquarters in Ocala, Florida, USA. If you are interested in business and/or advanced training, please call 352-401-0303.

We are very excited in presenting "The Traveler" (US Patent No. 5,871,463), the only portable colon hydrotherapy equipment in the world. Our patented, disposable speculums are available in 2 styles (Bulb and Straight) and also 2 sizes (Large and Small) size. They are manufactured and made of the finest quality medical grade plastic, which ensures patient comfort and safety.

We hope that you enjoy the packet and DVD on the Clearwater Colon Hydrotherapy systems and please feel free to contact us with any questions you may have.

*Your Health Matters,
Mary Ruth Baker, President -MA4612*



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON THERAPY EQUIPMENT & SUPPLIES INTRODUCING THE PPC-101 COLON THERAPY SYSTEM

FEATURES

PPC-101

- ◆ Moveable (on caster wheels)
- ◆ Quick disconnects for water, sewer & electrical line.
- ◆ All Stainless Steel Cabinet
- ◆ Size 41" high - 16" deep - 19" wide
- ◆ Weight 110 lbs

- ◆ Completely Closed System
- ◆ Sink & Faucet built in
- ◆ Drawer & Cabinet
- ◆ Self-Disinfectant System
- ◆ Natural Light Behind View Tube
- ◆ Pressure & Temperature Controlled with Safety Features
- ◆ Dependable and Safe Operation
- ◆ Specimen Collector Available



The PPC-101 is in compliance with all design and manufacturer guidelines set by Regulatory Agencies.

- ◆ Built -In Ultra Violet Light, with Water Filtration System with a 1 Micron Sediment and 0.5 Micron Carbon, 4 filters in the bottom cabinet.

Operational Educational Material are Available.




COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

The PPC-101 Clinical Model now comes with the Rust/Sediment, Carbon and Ultra-Violet water treatment system built-in!



Available: PPC-101 CLINICAL MODEL WITH SPECIMEN COLLECTOR





COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

PHYSICAL SPECIFICATIONS

Size: 41" Height
16" Depth
19" Width

Weight: 110 lbs.

Connections: Hot & Cold Water,
Electrical and Drainage Lines
Cabinet: Stainless Steel Construction
Wheels: Stainless Steel Casters

ACCESSORIES

- ◆ Disposable Speculums and Tubing
- ◆ Specimen Collector (Optional)



FEATURES

Clearwater Colon Hydrotherapy, Inc. offers a colonic irrigation system with a one year warranty.

The Professional Colonic Equipment is equipped with the following features:

- ◆ Completely Closed System
- ◆ Automatic Sanitation System
- ◆ Pressure & Temperature control & shut down valve
- ◆ Water Pressure Regulator - Pressure Relief System
- ◆ Storage Cabinet
- ◆ Two Quick Disconnects and lines for water hookup and sewage drainage lines and connections
- ◆ Natural Light Behind View Tube
- ◆ Built-In Rust/Sediment, Carbon & UV water treatment system

EXCLUSIVE FEATURES

- ◆ Drawer and cabinet for storage
- ◆ Sink & faucet for cleaning & disinfectant
- ◆ Portable Cabinet- located anywhere in the room
- ◆ Monitors Pressure throughout the treatment



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

All PPC's are user friendly. All models have the ability to be movable and portable.

WATER FILTERING SYSTEM

UV Water Filtration System (UV energy required for destruction of various organisms including bacteria and viruses)

- ◆ The Water Filtration System is specifically designed for Clearwater Colon Hydrotherapy systems.
- ◆ The four filter hot and cold unit is installed in the lower cabinet
- ◆ Standard features include UV lamp indicator with visual monitor
- ◆ Separate hot and cold filters
- ◆ The first stage of filtration is a 1 Micron Sediment filter, filtering sediment, rust, and dirt.
- ◆ The second stage of filtration is a 0.5 Micron Carbon, filtering taste and odor, lead reduction, cyst reduction (giardia, cryptosporidium).
- ◆ The third stage of filtration is a UV bulb, filtering bacteria and viruses.

PERFORMANCE SPECIFICATIONS

- ◆ Equipped with the latest filtering and UV water treatment systems
- ◆ Safe connect and disconnect of: Water, Sewer, and Electrical
- ◆ Portable
- ◆ Regulated Temperature and Pressure
- ◆ Accurate and sensitive indication of pressure and temperature
- ◆ Water in View Tube is full at all times for viewing and inspecting specimens
- ◆ Protection from contamination of the equipment and the client

INSTALLATION GUIDELINES

- ✓ Drain height of 12" or less
- ✓ Water pressure of 40 to 45
- ✓ Water heater of 40 gallons or greater
- ✓ Treatment table height 1 to 2 " higher than the view tube
- ✓ Hot and cold water supply lines
- ✓ 110/120 Vac. GFI electrical wall outlet



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON HYDROTHERAPY SYSTEM
PPC-101 (Professional Portable Colonic Equipment)
CLINICAL MODEL

MODEL PPC-101

INCLUDES: COLON HYDROTHERAPY SYSTEM PPC-101 WITH:

- Temperature and Pressure Safety Systems
- Self-Sanitizing Unit
- Built-In Sink & Faucet
- Manual & Installation, Operation (DVD)
- Built-In Rust/Sediment, Carbon & UV Water Treatment System

PPC-101 Clinical Model Without Specimen Collector: Call For Pricing

SUPPLY PACKAGE

(Required for the PPC-101 Equipment)

Supply Package Includes:

- Brushes (for cleaning view tubes)
- Disinfectant Solution
- Air Freshener
- Sanitation End Connectors (pack of twenty)
- Stainless Steel Cleaner
- Surface Disinfectant
- Two 1 Micron Rust/Sediment filter cartridges
- Two 0.5 Micron Carbon Block filter cartridges

Package Includes a 10% Discount:

\$ 134.64

Disposable Speculums needed - not included in supply package



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON HYDROTHERAPY SYSTEM
PPC-101 (Professional Portable Colonic Equipment)
CLINICAL MODEL WITH SPECIMEN COLLECTOR

MODEL PPC-101 WITH SPECIMEN COLLECTOR
INCLUDES: COLON HYDROTHERAPY SYSTEM PPC-101 WITH:

- Specimen Collector
- Temperature and Pressure Safety Systems
- Self-Sanitizing Unit
- Built-In Sink & Faucet
- Manual & Installation, Operation (DVD)
- Built-In Rust/Sediment, Carbon & UV Water Treatment System

PPC-101 Clinical Model With Specimen Collector: Call For Pricing

SUPPLY PACKAGE

(Required for the PPC-101 Equipment)

Supply Package Includes:

- Brushes (for cleaning view tubes)
- Disinfectant Solution
- Air Freshener
- Sanitation End Connectors (pack of twenty)
- Stainless Steel Cleaner
- Surface Disinfectant
- Two 1 Micron Rust/Sediment filter cartridges
- Two 0.5 Micron Carbon Block filter cartridges

Package Includes a 10% Discount:

\$ 134.64

Disposable Speculums needed - not included in supply package



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

PHYSICAL SPECIFICATIONS

Cabinet:	Carrying Case:
11" Height	25.3 Height
19" Depth	11.7 Depth
4.5"Width	18.6 Width
31.5 lbs.	21.5 lbs.

- Connections: Hot & Cold Water, Electrical and Drainage Lines
- Portable Carrying Case: On Wheels with a pull up handle for easy traveling.
- With Built-In Ultra-Violet System

OPTIONAL ACCESSORIES

- ◆ Disposable Speculums Kits
- ◆ Stainless Steel Cart with Sewage Tank
- ◆ Stand with Caster Wheels



FEATURES

Clearwater Colon Hydrotherapy, Inc. offers a colonic irrigation system with a one year warranty.

The Professional Portable Colonic Equipment is equipped with the following features:

- ◆ Completely Closed System
- ◆ Dual purpose capabilities: traveling or fixed locations
- ◆ Three options for Water Line connections
- ◆ Two options for Waste Line connections
- ◆ Pressure & Temperature control & shut down valve
- ◆ Water Pressure Regulator - Pressure Relief System
- ◆ Wall Bracket
- ◆ Built-In Rinse System after disinfecting
- ◆ Natural Light Behind View Tube
- ◆ Separate Water Filtration System
- ◆ UV Built into Equipment



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

TRAVELER DUAL PURPOSE CAPABILITIES

TRAVELING ABILITY

Water Line Connections (2 Options)

1. Single Sink Operations:
 - ◆ Single hot & cold water line combination. Line is affixed to the sink faucet for the hot and cold water supply.
 - ◆ The temperature is regulated at the sink faucet.
 - ◆ The colonic system mixing valve is then by passed, when using the combination hot and cold water supply inlet connection.
2. Dual sink operation:
 - ◆ Separate hot & cold water lines. Individual water lines are affixed to separate sink faucets. (If available)
 - ◆ One sink faucet is placed in full cold position and the other line is in full hot position and the lines are connected to the hot & cold inlet ports on the colonic equipment.
 - ◆ The temperature is regulated with the colonic equipment mixing valve.

PERFORMANCE SPECIFICATIONS

- ◆ Equipped with the latest filtration System
- ◆ UV built into Equipment
- ◆ Safe connect and disconnect of: Water, Sewer, and Electrical
- ◆ Portable
- ◆ Regulated Temperature and Pressure
- ◆ Accurate and sensitive indication of pressure and temperature
- ◆ Water in View Tube is full at all times for viewing
- ◆ Protection from contamination to the equipment and the client

INSTALLATION GUIDELINES

- ✓ Drain height of 12" or less
- ✓ Water pressure of 40 to 45
- ✓ Water heater of 40 gallons or greater (greater is better)
- ✓ Treatment table height 1 to 2 " higher than the view tube
- ✓ Hot and cold water supply lines
- ✓ 110/120 Vac. GFI electrical wall outlet



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON HYDROTHERAPY SYSTEM
PPC-101 (Professional Portable Colonic Equipment)
THE TRAVELER

MODEL PPC-101 "THE TRAVELER"

INCLUDES: COLON HYDROTHERAPY SYSTEM PPC-101 WITH:

- Temperature and Pressure Safety Systems
- Self-Sanitizing Unit
- Carrying Case on wheels with pull up handle
- Manual & Installation/Operation (DVD)
- Built-in Ultra-Violet System, including a separate unit with Hot & Cold - 4 filter stainless steel Water Filtration Treatment System

PPC-101 "The Traveler" Model (U.S. Patent No. 5,871,463): Call Us

SUPPLY PACKAGE

(Required for the PPC-101 Equipment)

Supply Package Includes:

- Brushes (for cleaning view tubes)
- Disinfectant Solution
- Air Freshener
- Sanitation End Connectors (pack of twenty)
- Stainless Steel Cleaner
- Surface Disinfectant
- Two 1 Micron Rust/Sediment filter cartridges
- Two 0.5 Micron Carbon Block filter cartridges

Package Includes a 10% Discount:

\$ 134.64

Disposable Speculums needed - not included in supply package



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON HYDROTHERAPY SYSTEM
PPC-101 (Professional Portable Colonic Equipment)
“THE TRAVELER” With Options
Sewage Tank Stand and
Stainless Steel Stand

**Stainless Steel Cart with Sewage Tank
for Traveler**

- ◆ 33 Gallon Septic Tank - Drainage System
- ◆ with Handicap Toilet Height
- ◆ Lower Storage Tray
- ◆ Adjustable Height for Traveler (Range 31.5 to 40.0 inches)
- ◆ On Locking Caster Wheels (Wt. 113 lbs)



Stainless Steel Stand for Traveler

- ◆ Meets Required Height of 28"
- ◆ Lower Storage Tray for Filters
- ◆ On Locking Caster Wheels





COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Water Filtration System Supplies

1. The first stage of filtration is a 1 Micron Sediment filter, filtering sediment, rust, and dirt.
2. The second stage of filtration is a 0.5 Micron Carbon, filtering taste and odor, lead reduction, cyst reduction (giardia, cryptosporidium).
3. The third stage of filtration is a UV bulb, filtering bacteria and viruses.

Granular Activated Carbon Cartridges

- ◆ 5 inch and 10 inch



Fiber Sediment Cartridges

- ◆ 5 inch and 10 inch



Carbon Filter (Gac Inline) for Other Manufacturer Units

- ◆ 10 inch





COLON HYDROTHERAPY, INC.®


3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

NEW SPECULUM INFLOW WATER LINE CONNECT AND DISCONNECT DIRECTIONS - D/C

- 1) On the left side of the unit you will find a male brass style fitting with a white water line disconnect facing outward - as shown in the photo.



- 2) **To Connect:** Take the white disconnect and connect it to the line by a simple alignment then push and twist clockwise (to the right). You should hear and feel it click into place.



- 3) **To Disconnect the inflow line:** Twist the white disconnect counter clockwise (to the left) with a 1/4 turn and then pull straight out.



PLEASE NOTE:

This quick and easy connector device will be on each and every Speculum kit line. Be sure when ordering to indicate that you have the new DC connector.

FOR OLDER CLEARWATER UNITS: You can convert your equipment to the new D/C connector. Call us for the new conversion kit, which will be shipped free of charge with your next order. It will save ware and tare on your unit and you will no longer have to cut off the hose with a sharp object.



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Disposable Speculums

Available in Large and Small sizes, and Bulb and Straight styles



FEATURES

Large Straight or **Small Straight**
Large Bulb or **Small Bulb**
48", 42" or 36" Waste Hose length
1/4" or 3/8" Inflow Hose
Standard Hose or G-Hose

Designed for easier insertion and greater client comfort. Obturator is also designed with a longer thumb tab for ease in handling.
Lubricant packet (net wt. 5g) included.

BIO-COMPATIBILITY



TEST RESULTS:

An independent laboratory study was conducted to determine bio-compatibility on the speculums and hosing. The laboratory was looking for harmful: chemicals and toxicity that would be leached from the plastics itself. There was no findings of resins or toxicity found that could be derived by contact from this material. It was given a Superior Medical Rating.

BULB STYLE

(PATENT NO. US D452004S)

STRAIGHT STYLE





COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

NEVER RUN OUT OF STOCK WITH OUR AUTOSHIP PROGRAM

According to several management studies, the single most important characteristic of an effective office is the ability to manage time and worries.

Ask yourself these questions:

- ☹ Are you running out of stock?
- ☹ Do you have a hard time meeting deadlines?
- ☹ Do you hate keeping up with your inventory?
- ☹ Do you call to overnight items, because you're running short?
- ☹ Do you worry about canceling appointments, because you do not have what you need?
- ☹ Do you hastily call to order stock, pressed for time, only to remember something more after you already hung up the phone?
- ☹ Are you busy with patients and do not have the time to call in the orders?
- ☹ Were you ever "out of business" because your stock hasn't arrived or you forgot to place that order?
- ☹ Do you check with everyone to make sure someone has placed your order?

If you answered YES to even one or more of these questions, you're sure to benefit from our popular, effective method to ship what you need, and have it there, when you need it.

The program I just introduced is our **AUTO-SHIP Program**.

- ☺ You choose what you want shipped.
- ☺ When you want it shipped.
- ☺ How much you want shipped.



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

But the best benefit is☺.... you get a discount for having us help you out.

You now have peace of mind, NO stress, NEVER run out. Results: SAVE TIME and MONEY.

Use the enclosed “AutoShip Program Agreement” in the pocket of this catalog to tailor your Auto-Shipment needs and drop it in the mail, fax it, or give us a call and we will help you get started.

Some suggestions along with your speculum order are, ordering a case of Underpads every other order, or adding Filters and gloves as needed. We can also take a look at your ordering history, which can help you in your decision.

We want to show you how to eliminate an unnecessary interruptions, so we suggest that you save time by making time to plan your Auto-Shipment right now, while you still have a mind to save time. Your business will be exploding!!

☺ We Have It All Under Control.





COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Therapy Table & Accessories

Adjustable Electric Therapy Table



ADVANTAGES:

(Qualifies for Wheel Chair Accessible Tax Credit)

- Lowers to Wheel Chair Height
- Adjustment for Patient and Therapist Body Mechanics
- A Fully Equipped Bodywork Table
- Minimal Maintenance

FEATURES:

- Available in 30" wide x 6' long
- 10" Extension can be added to the table
- Two Designs (Flat or Tilt Back)
- Upholstery color is Buff
- Frame color is a textured ice white epoxy paint

SPECIFICATIONS:

- Top of the pad height at lowest point is 17"
- Top of the pad height at highest point is 30"
- Pad is 3" of med density foam
- Air operated foot controls
- Standard 115 VAC, 4 AMPS
- Standard 8 ft power cord



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Light Weight Portable Table and Accessories



- Upholstery Standard Color - Blue or Black (In Stock)
- Other Upholstery Colors - By Special Order Only
- **Table Cart Carrier** - Weight 9 lbs, but will hold 180 lbs
- 27" x 14" x 7" When Closed
- Enjoy the ease of rolling your table where ever you go.



- **Pillows: Provide the correct level of support for prone, supine and side-lying positions**
 - ▶ **Wedge** - With Matching Vinyl Cover
 - ▶ **Round Bolster** - 6" x 6" x 26"
 - ▶ **Large Round Bolster** - 9" x 9" x 26"
(Support under the knees, ankles or use for side posture support)





COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Therapist's Chair and Accessories

- Designed for the professional therapists
- Designed for comfort of sitting for long periods
- Seat in your choice of color
- Easy rolling on caster wheels



- Lightweight, yet sturdy
- "No-slip, no-tip" design
- Durable top, covered with black ribbed rubber and trimmed with grey molding



- Waterproof Table Cover
- Protects your table from moisture and grime
- No liquids seep into your table's foam



- Table Warmer with Thermostat Controls
- 2 styles - Basic and Deluxe
- Keep clients warm, comfortable, and happy during treatment





COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474


Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com


Colon Therapy Charts



COLON MODEL



COLON THERAPY CHART



THE DIGESTIVE SYSTEM



BLUEPRINT FOR HEALTH




COLON HYDROTHERAPY, INC.®


3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191


www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com




DISEASES OF DIGESTIVE SYSTEM



GASTROESOPHAGEAL DISORDERS



BOWEL REFLEX POINTS



(IBD) INFLAMMATORY BOWEL DIS




COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com



UNDERSTANDING IRRITABLE BOWEL SYNDROME (IBS)



STAINLESS STEEL VIEW TUBE BRUSH HOLDER

GOWNS

- ◆ Custom Made
- ◆ Small, Medium, Large, X-Large, and XX-Large
- ◆ Two different Lengths (To the Calves, or ankles)
- ◆ Fade-Resistant, Wrinkle-Free
- ◆ Cotton/Poly Seersucker Material
- ◆ Soft Elegant Pastel Stripes
- ◆ Velcro Closures







COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com



THESE TWO GREAT NON-GLUTEN RECIPE BOOKS ARE A MUST FOR YOUR COOKBOOK LIBRARY! BY GWYNN PALMER



SUPPLY PACKET



SPECULUM SAMPLE PACKET



EXAM GLOVES



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Attention Equipment and/or Speculum Buyers:

In accordance with Federal Regulations, it is a requirement that a letter be placed in the file of the equipment and/or purchasing speculums. When you purchase our equipment, we will need one for the equipment and one for the ordering of the disposable speculums. A sample letter has been enclosed for you to follow exactly with no variations for both in one script.

“Caution: Federal law restricts this device to sale by or on the order of a physician.” 21CFR801.109

This regulation applies to all Class II Colon Therapy Equipment and Disposable Speculums by every Manufacturer both United States and Foreign Countries.

We also have a dual device listing on our disposables speculums for Class I devices, which has different package labeling. To ensure you are in compliance with the Federal Law, we will need to be advised of the class of device that you are using so that we may provide you with the proper package labeling.

Sample Letter

Letter Head with Doctor’s Name, Address and Phone number

Clearwater Colon Hydrotherapy, Inc.
3145 S.W. 74th Terrace
Ocala, Florida 34474

Subject: Order to Purchase
Re: PPC-101 Colon Hydrotherapy Unit and/or Disposable Speculums

In Accordance with the Federal Code of Regulations 21 Food and Drugs, Parts 800-1299 as stated below.

“Caution: Federal law restricts this device to sale by or on the order of a physician.”
21CFR801.109

This letter is to serve as a order for the following:

Person: _____ Clinic: _____ to
purchase PPC-101 Colon Hydrotherapy Unit and/or Disposable Speculums/Kits.

Sincerely, **Doctor’s Signature**
Printed Name and License #



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Boost your Business by Using Video E-Mail

Did you know....that Over 40 billion E-mails are sent daily, yet less than 1% of internet users have experienced Video E-mail.

DON'T TYPE TALK

DON'T EXPLAIN HOW TO DO SOMETHING, SHOW THEM LIVE

This is what we use at Clearwater Colon Hydrotherapy. You probably have already opened one or more of our promotional and educational videos. Individuals and businesses alike can take advantage of the Video Communication Solution in promoting your Business.

Everyone can easily stay in touch with clients and customers as well as family and friends around the world, sharing videos and personal thoughts and messages, in a way no other means of communication can offer. No more boring black and white text E-mails. Promote your business with ease, no experience necessary. Create personal and business E-mails, and if you don't like the camera, use their templates. Imagine never missing another incredible moment, no matter where you are. It is also one of the best unsurpassed marketing tool for promoting your business or service, anywhere in the world. You can keep in contact with existing customers, introduce yourself to new clients, send out Video Newsletters, let people know about special offers, send thank you's, remind them of their upcoming appointment,... the list of potential uses goes on and on.

You can even deliver LIVE presentations anytime, anywhere, to an unlimited audience. You can conduct your meetings online easily and efficiently. Even set up your Webinars easily in minutes to unlimited viewers worldwide!

This program will change the way you communicate your business.

The history of communication shows ... we started with RADIO - the spoken word then along came TELEVISION, which added pictures to the spoken word

Then, E-Mail, which is the written word

and now VIDEO E-MAIL, which simply is adding Video to E-mail.

If you would like additional information, on how to get started, Clearwater can send you a couple of "Product Videos" for you to view, like showing you how you can create the Talking Video - "Promoting Your Business In Just 60 seconds". It's so simple and easy, yet you have a Management Team that has over 150 years of combined experience in Marketing, Management and Direct Selling, using all their Internet Technology behind you. For More Information:

Please give us a call at (352) 401-0303 or E-Mail us at sales@colonhydrotherapy.com



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Product Information and Price List

ULTRA-VIOLET WATER TREATMENT SYSTEMS

Ultra-Violet Lamp (Replacement Bulbs 10000 Hour UV Bulb \$49.95

PROFESSIONAL HYDROTHERAPY MODELS

PPC-101 The Traveler Call For Pricing
 PPC-101 Clinical Call For Pricing
 PPC-101 Clinical with Specimen Collector Call For Pricing

Options For Traveler:

Stainless Steel Cart with Sewage Tank Call For Pricing
 Stainless Steel Mobile Table Call For Pricing

DISPOSABLE SPECULUMS

Large & Small BULB Styles \$3.95
 Large & Small STRAIGHT Styles \$3.95
 Kits with D/C Connectors \$4.45
 Kits Through the Auto-Shipment Program \$3.65
 Speculum Sample Pack - Standard \$14.60
 Speculum Sample Pack - D/C \$16.60
 (The Sample Packets Consists of: 4 Complete Kits - 1 Kit for Each Size and Style of
 Speculums)

WATER FILTER REPLACEMENT CARTRIDGES

Granular Activated Carbon Cartridges (0.5 Micron Taste/Order/Chemical)
 10 Inch \$19.95
 5 Inch \$19.50
 Fiber Sediment Cartridges (Rust/Sediment 1 Micron)
 10 Inch \$6.00
 5 Inch \$5.00
 Carbon Filter (GAC Inline) For Other Manufacturer Units
 10 Inch \$9.95



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

CLINICAL & PORTABLE COLON THERAPY TABLES AND ACCESSORIES

Light Weight Portable Table (27 Pounds)	\$459.00
Adjustable Electric Therapy Table Flat Top	\$1,399.00
Adjustable Electric With Manual Tilt Option	\$1,499.00
Table Cart Carrier	\$89.00
Plastic Sheet (Table Cover)	\$15.95
Table Cover Warmer (Basic)	\$69.00
Table Cover Warmer (Deluxe)	\$99.00
Therapist's Chair	\$425.00
Pillow Wedge (with matching cover)	\$130.00
Bolster Full Round (with matching cover)	\$45.95
Bolster Jumbo Round (with matching cover)	\$63.00
Step Stool (Medical Grade Construction)	\$59.00
Toilet Step	\$59.00

COLON THERAPY SUPPLIES

Supply Packet (10% Discount Applied)	\$134.64
(Packet Consists of: - Pack of two View Tube Brushes, Disinfectant Solution, Order Neutralizer/Cleaner, Pack of 20 Sanitation End Connectors, Stainless Steel Cleaner, Surface Disinfectant, 2 Sediment Filters, and 2 Carbon Filters)	
Surface Disinfectant Decontaminate Cleaner	\$12.50
Disinfectant Solution (Approved for Endoscopic Equipment 1 Gallon)	\$28.55
Sanitation End Connectors (Pack of 20)	\$3.25
Stainless Steel Cleaner and Polish (1 Spray Can 16 oz)	\$17.35
Odor Neutralizer/Cleaner Spray (16 oz)	\$14.95
View Tube Brushes (Pack of two 24" Length)	\$24.00
Stainless Steel Brush Holder	\$45.00
View Tube for PPC-101 Clinical Model (Quantity 1)	\$50.00
View Tube for PPC-101 Traveler Model (Quantity 1)	\$50.00
Light Bulb for View Tube Assembly (Quantity 1) 9000 hr Natural Lighting	\$18.00



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Exam Gloves - Powder Free, Nitrile Gloves	\$8.85
(Small, Medium, Large and X-Large)	
Underpads (Large Case of 150)	\$49.98
Disposable Footies	\$39.00
Disposable Gowns (Case of 50)	\$26.55
Gowns (See Sizes Below)	
Large - 42"	\$18.95
X-Large - 42"	\$20.95
XX-Large - 42"	\$22.95
Large - (Custom 51" Size to Ankle)	\$24.95
Extra Large - (Custom 51" Size to Ankle)	\$26.95
Lubricant (Tube)	\$3.75
Hemorrhoid Relief Cream	\$16.95

*Recommended Reading for Certification. These are Excellent References to have in your library. Items Marked with an * are Suggested.*

BOOKS

* Office Safety Manual (Set to OSHA Standards for Colon Therapy)	
Author: Mary Ruth Baker	\$35.00
* Dr. Protocol Manual - Author: Mary Ruth Baker	\$35.00
* Staying Healthy with the Seasons	\$15.95
* The CIBA Collection of Medical Illustrations	
Upper Digestive Tract	\$105.00
Lower Digestive Tract	\$105.00
Eat Well and Live (Gluten Free - Author: Gwynn Palmer)	\$21.95
Having It All (Gluten Free - Author: Gwynn Palmer)	\$22.95
Colon Health: The Key to a Vibrant Life	\$9.95
8 Weeks to Women's Wellness	\$16.95



COLON HYDROTHERAPY, INC.®

3145 SW 74th Terrace, Ocala, FL 34474

Ph: 352-401-0303 ~ Fax: 352-401-9197 ~ Toll Free Order Line: 888-869-6191

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

CD's

Medicine Woman (CD) Part I	\$15.95
Medicine Woman (CD) Part II	\$15.95
Medicine Woman (CD) Part III	\$15.95
Wheat Belly - (Set of 6 CD's)	\$29.95
Why People Don't Heal and How They Can - (Set of 6 CD's)	\$24.95

MODELS

Colon Model	\$82.80
-------------------	---------

CHARTS

Colon Therapy (Chart of the Colon)	\$19.95
Digestive System	\$19.95
Blueprint for Health Chart	\$19.95
Diseases of the Digestive System	\$19.95
Gastroesophageal Disorders	\$19.95
Bowel Reflex Points	\$14.95
Inflammatory Bowel Disease (IBD)	\$19.95
Understanding Irritable Bowel Syndrome (IBS)	\$19.95

Quick Check List - Morning Set-up, Therapy & Disinfecting the Unit After Use	
For Clinical	\$3.95
For Traveler	\$3.95

PLEASE NOTE: Prices are subject to change without notice and do not include shipping charges.