

3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com



2023 CATALOG

PRESS RELEASE - 1st In the Industry

Cytotoxicity Tested and Titanium Oxide Free Speculums and DEHP Free Tubing California Proposition 65 Compliant

CLEARWATER IS GOING




GREEN

Clearwater Colon Hydrotherapy, Inc. Product Catalog - Rev. 9.3 3-29-2023



3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com



LABELING

ADEQUATE DIRECTIONS FOR USE:

Under no circumstances should the patient be left unattended. 21CFR801.109

INDICATION FOR USE STATEMENT:

The indication for use of this device must be restricted to colon cleansing when medically indicated, such as before radiological or endoscopic examination.

21CFR876.5220

PRESCRIPTION STATEMENT:

"Caution: Federal law restricts this device to sale by or on the order of a physician."

21CFR801.109



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON PREP EQUIPMENT AND SUPPLIES

TABLE OF CONTENTS

Clearwater Colon Hydrotherapy Equipment	
Labeling	1
Benefits of Our Equipments	5
Background	
PPC-101 Clinical Model with and w/o Specimen Collector	7
PPC-101 Clinical Model with Built In UV Filtration System	8
PPC-101 Clinical Features	9
PPC-101 "The Traveler" Features	
PPC-101 "The Traveler" Dual Purpose	14
PPC-101 "The Traveler" with Sewage Tank Stand	
"The Traveler" Stainless Steel Stand	16
Water Filtration System	
Granular Activated Carbon Cartridges (10 Inch)	17
(1 Micron Taste/Odor/Chemical) (5 Inch)	17
Fiber Cartridges (10 Inch)	17
(Rust/Sediment 1 Micron (5 Inch)	
Carbon Filter (Gac Inline) for Other Manufacturer Units (10 Inch)	17
Disposable Speculums	
Speculum Inflow water Line Disconnects	18
Bio-Compatibility Test Results	19
	19
AutoShip Program	20



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Adjustable Electric Therapy Table	22
Light Weight Portable Tables	23
Fable Cart Carrier	23
Pillow Wedge with Vinyl Cover	23
Bolster Round Supports (Knees/Ankles)	23
Therapist's Chair	24
±	
Step Stool	
Γable Cover	24
Fable Cover Warmer	24

Colon Therapy Charts

Therapy Table & Accessories

Colon Therapy	25
Digestive System	
Blueprint for Health Chart	
Diseases of the Digestive System	
Gastroesophageal Disorders and Digestive Anatomy	
Bowel Reflex Points	26
Inflammatory Bowel Disease (IBD)	26
Understanding Irritable Bowel Syndrome (IBS)	



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Products & Supplies	
Stainless Steel View Tube Brush Holder	
Books (Complete List in Product Info & Price Section pg 33)28Supply Packet28Speculum Sample Pack28Gloves28	8
Federal Code of Regulations on Purchasing Device and Speculums Sample Letter/Script on Doctors Letterhead	•
Boosting Your Business - Using Video E-Mail Video E-Mail Program)
Product Information and Price List Item Price List	1
Office Order Form Fax Order Form for Purchases	5



SOME OF THE REASONS WHY OUR EQUIPMENT WILL GREATLY BENEFIT YOU AND YOUR PATIENTS!

Equipment used for cleansing the colon - such as, before radiological or endoscopic examinations. Used for preparation for a colonoscopy.

- ♦ Designed for use in Gastroenterology Clinic, Hospitals & Nursing Homes
- **♦** Cuts down on bad preps
- ♦ Portable can go to bedside of Colonoscopy table
- ♦ Sewage Tank Stand is great when plumbing can not be done

Abstract on Colon Hydrotherapy for Pre-endoscopy Preparation (Copy of Abstract upon request)

Do you have bad preps?

- ♦ Our Equipment monitors the Colon pressure throughout the whole session
- **♦** Completely closed system
- **♦** All Stainless Steel Cabinets
- **♦** Pressure & Temperature Controlled with safety features
- ♦ Water in View Tube is full at all times for viewing and inspecting specimens
- ♦ Safe connects and disconnects of Water, Sewer, and Electrical
- **♦** Options for Water Line and Waste Line connections
- ♦ Quick and easy connector device on each and every Speculum kit line
- ♦ Disposable Speculums that comes in different sizes and styles
- ♦ Speculums are designed for easier insertion and greater client comfort
- ♦ Bio-Compatibility test results study given Superior Medical Rating
- ♦ Free 2 full days of Training with purchase of a unit

Do you have different locations?

- ♦ Our Equipment is the answer. We have the Patent on Portability
- **♦** Patent on the ONLY Speculum (rectal nozzles) that are medical grade plastic and designed anatomically correct.
- ♦ ONLY inflow water line that is DEHP FREE (Important for your sensitive patients)



Clearwater Colon Hydrotherapy, was established in 1992. Mary Ruth Baker, owner, has been in the field of Colon Hydrotherapy for over 35 years. Her commitment has always been to have the highest standard of quality and client care. These high standards are now available in Clearwater Colon Hydrotherapy systems. These systems have the best mechanical and operational design, which makes them stand alone as the most efficient colon prep for Colonoscopy, colon cleansing equipment that is on the market today.

Mary Ruth is anxious to teach and support you in the operations of the Clearwater Colon Hydrotherapy equipment and also with the knowledge of the business operations towards setting up a complaint office.

An installation manual & video is provided when you purchase a unit, that demonstrates the complete installation process. A two day training seminar on the installation, operation and maintenance of the equipment is provided free of charge, if you schedule your training within 30 days of your initial equipment purchase, at our corporate headquarters in Ocala, Florida, USA. If you are interested in business and/or advanced training, please call 352-401-0303.

We are very excited in presenting "The Traveler" (US Patent No. 5,871,463), the only portable colon hydrotherapy equipment in the world. Our patented, disposable speculums are available in 2 styles (Bulb and Straight) and also 2 sizes (Large and Small)size. They are manufactured and made of the finest quality medical grade plastic, which ensures patient comfort and safety.

We hope that you enjoy the packet and Video on the Clearwater Colon Hydrotherapy systems and please feel free to contact us with any questions you may have.

Your Health Matters, Mary Ruth Baker, President - MA4612



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON THERAPY EQUIPMENT & SUPPLIES INTRODUCING THE PPC-101 COLON THERAPY SYSTEM

FEATURES

PPC-101

- Moveable (on caster wheels)
- Quick disconnects for water, sewer & electrical line.
- All Stainless Steel Cabinet
- Size 41" high 16" deep 19" wide
- Weight 110 lbs
- Completely Closed System
- Sink & Faucet built in
- Drawer & Cabinet
- Self-Disinfectant System
- Natural Light Behind View Tube
- Pressure & Temperature Controlled with Safety Features
- Dependable and Safe Operation
- Specimen Collector Available



The PPC-101 is in compliance with all design and manufacturer guidelines set by Regulatory Agencies.

Built -In Ultra Violet Light, with Water Filtration System with a 1 Micron Sediment and 0.5 Micron Carbon, 4 filters in the bottom cabinet.

Operational Educational Material are Available.



The PPC-101 Clinical Model now comes with the Rust/Sediment, Carbon and Ultra-Violet water treatment system built-in!



Ultra-Violet treatment system

Rust/Sediment Filters

Activated carbon filters (removes taste/odor/chemicals)

Available: PPC-101 CLINICAL MODEL WITH SPECIMEN COLLECTOR







3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

PHYSICAL SPECIFICATIONS

Size: 41" Height

16" Depth 19"Width

Weight: 110 lbs.

Connections: Hot & Cold Water, Electrical and Drainage Lines

Cabinet: Stainless Steel Construction

Wheels: Stainless Steel Casters

ACCESSORIES

- Disposable Speculums and Tubing
- ♦ Specimen Collector (Optional)



FEATURES

Clearwater Colon Hydrotherapy, Inc. offers a colonic irrigation system with a one year warranty.

The Professional Colonic Equipment is equipped with the following features:

- ♦ Completely Closed System
- ♦ Automatic Sanitation System
- ♦ Pressure & Temperature control & shut down valve
- Water Pressure Regulator -Pressure Relief System
- ♦ Storage Cabinet
- ♦ Two Quick Disconnects and lines for water hookup and sewage drainage lines and connections
- Natural Light Behind View Tube
- ♦ Built-In Rust/Sediment, Carbon & UV water treatment system

EXCLUSIVE FEATURES

- Drawer and cabinet for storage
- Sink & faucet for cleaning & disinfectant
- ♦ Portable Cabinet- located anywhere in the room
- ♦ Monitors Pressure throughout the treatment



All PPC's are user friendly. All models have the ability to be movable and portable.

WATER FILTERING SYSTEM

UV Water Filtration System (UV energy required for destruction of various organisms including bacteria and viruses)

- ♦ The Water Filtration System is specifically designed for Clearwater Colon Hydrotherapy systems.
- ♦ The four filter hot and cold unit is installed in the lower cabinet
- ♦ Standard features include UV lamp indicator with visual monitor
- ♦ Separate hot and cold filters
- ♦ The first stage of filtration is a 1 Micron Sediment filter, filtering sediment, rust, and dirt.
- ♦ The second stage of filtration is a 0.5 Micron Carbon, filtering taste and odor, lead reduction, cyst reduction (giardia, crytosporidium).
- ♦ The third stage of filtration is a UV bulb, filtering bacteria and viruses.

PERFORMANCE SPECIFICATIONS

- ♦ Equipped with the latest filtering and UV water treatment systems
- ♦ Safe connect and disconnect of: Water, Sewer, and Electrical
- **♦** Portable
- Regulated Temperature and Pressure
- ♦ Accurate and sensitive indication of pressure and temperature
- ♦ Water in View Tube is full at all times for viewing and inspecting specimens
- ♦ Protection from contamination of the equipment and the client

INSTALLATION GUIDELINES

- ✓ Drain height of 12" or less
- ✓ Water pressure of 40 to 45
- ✓ Water heater of 40 gallons or greater
- ✓ Treatment table height 1 to 2" higher than the view tube
- ✓ Hot and cold water supply lines
- ✓ 110/120 standard wall outlet



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON HYDROTHERAPY SYSTEM PPC-101 (Professional Portable Colonic Equipment) CLINICAL MODEL

MODEL PPC-101 INCLUDES: COLON HYDROTHERAPY SYSTEM PPC-101 WITH:

- Temperature and Pressure Safety Systems
- Self-Sanitizing Unit
- Built-In Sink & Faucet
- Manual & Installation & Video
- Built-In Rust/Sediment, and Carbon & UV Water Treatment System

PPC-101 Clinical Model Without Specimen Collector: Call In

SUPPLY PACKAGE

(Required for the PPC-101 Equipment)

Supply Package Includes:

- Brushes (for cleaning view tubes)
- **Disinfectant Solution**
- Air Freshener
- Sanitation End Connectors (pack of twenty)
- Stainless Steel Cleaner
- Surface Disinfectant
- Two 1 Micron Rust/Sediment filter cartridges
- Two 0.5 Micron Carbon Block filter cartridges

Package Includes a 10% Discount:

\$ 155.64

Disposable Speculums needed - not included in supply package



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON HYDROTHERAPY SYSTEM PPC-101 (Professional Portable Colonic Equipment) CLINICAL MODEL WITH SPECIMEN COLLECTOR

MODEL PPC-101 WITH SPECIMEN COLLECTOR INCLUDES: COLON HYDROTHERAPY SYSTEM PPC-101 WITH:

- Specimen Collector
- Temperature and Pressure Safety Systems
- Self-Sanitizing Unit
- Built-In Sink & Faucet
- Manual & Installation & Video
- Built-In Rust/Sediment, Carbon & UV Water **Treatment System**

PPC-101 Clinical Model With Specimen Collector: Call In

SUPPLY PACKAGE

(Required for the PPC-101 Equipment)

Supply Package Includes:

- Brushes (for cleaning view tubes)
- **Disinfectant Solution**
- Air Freshener
- Sanitation End Connectors (pack of twenty)
- Stainless Steel Cleaner
- Surface Disinfectant
- Two 1 Micron Rust/Sediment filter cartridges
- Two 0.5 Micron Carbon Block filter cartridges

Package Includes a 10% Discount:

\$ 155.64



3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Disposable Speculums needed - not included in supply package

PHYSICAL SPECIFICATIONS

Cabinet: Carrying Case: 11" Height 25.3 Height 4.5" Depth 11.7 Depth 18.6 Width 31.5 lbs. 21.5 lbs.

- Connections: Hot & Cold Water, Electrical and Drainage Lines
- Portable Carrying Case: On Wheels with a pull up handle for easy traveling.
- With Built-In Ultra-Violet System

OPTIONAL ACCESSORIES

- **♦**Disposable Speculums Kits
 - ♦ Stainless Steel Cart with Sewage Tank



♦ Stand with Caster Wheels **FEATURES**

Clearwater Colon Hydrotherapy, Inc. offers a colonic irrigation system with a one year warranty.

The Professional Portable Colonic Equipment is equipped with the following features:

- ♦ Completely Closed System
- Dual purpose capabilities: traveling or fixed locations
- ♦ Three options for Water Line connections
- ♦ Two options for Waste Line connections
- ♦ Pressure & Temperature control & shut down valve
- ♦ Water Pressure Regulator Pressure Relief System
- ♦ Wall Bracket
- ♦ Built-In Rinse System after disinfecting
- ♦ Natural Light Behind View Tube
- ♦ Separate Water Filtration System



3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

♦ UV Built into Equipment

TRAVELER DUAL PURPOSE CAPABILITIES

TRAVELING ABILITY

Water Line Connections (2 Options)

- 1. Single Sink Operations:
- ♦ Single hot & cold water line combination. Line is affixed to the sink faucet for the hot and cold water supply.
- ♦ The temperature is regulated at the sink faucet.
- ♦ The colonic system mixing valve is then by passed, when using the combination hot and cold water supply inlet connection.
- 2. Dual sink operation:
- ♦ Separate hot & cold water lines. Individual water lines are affixed to separate sink faucets. (If available)
- ♦ One sink faucet is placed in full cold position and the other line is in full hot position and the lines are connected to the hot & cold inlet ports on the colonic equipment.
- ♦ The temperature is regulated

with the colonic equipment mixing valve.

PERFORMANCE SPECIFICATIONS

- ♦ Equipped with the latest filtration System
- ♦ UV built into Equipment
- ♦ Safe connect and disconnect of Water, Sewer, and Electrical
- **♦** Portable
- Regulated Temperature and Pressure
- ♦ Accurate and sensitive indication of pressure and temperature
- ♦ Water in View Tube is full at all times for viewing
- ♦ Protection from contamination to the equipment and the client

INSTALLATION GUIDELINES

- ✓ Drain height of 12" or less
- ✓ Water pressure of 40 to 45
- ✓ Water heater of 40 gallons or greater (greater is better)
- ✓ Treatment table height 1 to 2 " higher than the view tube
- ✓ Hot and cold water supply lines
- ✓ 110/120 Volt electrical wall outlet



COLON HYDROTHERAPY SYSTEM PPC-101 (Professional Portable Colonic Equipment) THE TRAVELER

MODEL PPC-101 "THE TRAVELER" INCLUDES: COLON HYDROTHERAPY SYSTEM PPC-101 WITH:

- Temperature and Pressure Safety Systems
- Self-Sanitizing Unit
- Carrying Case on wheels with pull up handle
- Manual & Installation/Operation & Video
- Built-in Ultra-Violet System, including a separate unit with Hot & Cold - 4 filter stainless steel Water Filtration Treatment System

PPC-101 "The Traveler" Model (U.S. Patent No. 5,871,463): Call In

SUPPLY PACKAGE

(Required for the PPC-101 Equipment)

Supply Package Includes:

- Brushes (for cleaning view tubes)
- Disinfectant Solution
- Air Freshener
- Sanitation End Connectors (pack of twenty)
- Stainless Steel Cleaner
- Surface Disinfectant
- Two 1 Micron Rust/Sediment filter cartridges
- Two 0.5 Micron Carbon Block filter cartridges

Package Includes a 10% Discount:

\$ 155.64

Disposable Speculums needed - not included in supply package

Clearwater Colon Hydrotherapy, Inc. Product Catalog - Rev. 9.3 3-29-2023



3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

COLON HYDROTHERAPY SYSTEM PPC-101 (Professional Portable Colonic Equipment) "THE TRAVELER" With Options

Sewage Tank Stand and Stainless Steel Stand

Stainless Steel Cart with Sewage Tank for Traveler

- ♦ 33 Gallon Septic Tank Drainage System
- **♦** with Handicap Toilet Height
- **♦** Lower Storage Tray
- ♦ Adjustable Height for Traveler (Range 31.5 to 40.0 inches)
- ♦ On Locking Caster Wheels (Wt.113 lbs)

Stainless Steel Stand for Traveler

- ♦ Meets Required Height of 28"
- **♦** Lower Storage Tray for Filters
- **♦** On Locking Caster Wheels







3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Water Filtration System Supplies

- 1. The first stage of filtration is a 1 Micron Sediment filter, filtering sediment, rust, and dirt.
- 2. The second stage of filtration is a 0.5 Micron Carbon, filtering taste and odor, lead reduction, cyst reduction (giardia, crytosporidium).
- 3. The third stage of filtration is a UV bulb, filtering bacteria and viruses.

Granular Activated Carbon Cartridges ♦ 5 inch and 10 inch



Fiber Sediment Cartridges

♦ 5 inch and 10 inch



Carbon Filter (Gac Inline) for Other Manufacturer Units

♦ 10 inch





3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

NEW SPECULUM INFLOW WATER LINE CONNECT AND DISCONNECT DIRECTIONS - D/C

1) On the left side of the unit you will find a male brass style fitting Female disconnect facing outward - as shown in the photo.



2) **To Connect:** Take the white disconnect and connect it by a simply alignment then push in. You should hear and feel it click into place.



3) **To Disconnect the inflow line:** Push on the Female disconnect inward and the fitting will pop out.



PLEASE NOTE:

This quick and easy connector device will be on each and every Speculum kit line. Be sure when ordering to indicate that you have the new DC connector. This more sanitary feature for the equipment.

FOR ANY CLEARWATER UNITS: You can convert your equipment to the new CCH D/C connector. Call us and indicate that you want the new connector, fits all units ever built. It will save ware and tare on your unit and you will no longer have to cut off the hose with a sharp object.



3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Disposable Speculums

Available in Large, Medi and Small sizes, Bulb and Straight styles





FEATURES

BIO-COMPATIBILITY TEST RESULTS:

Large Straight or Small Straight, Medi
Straight, Large Bulb or Small Bulb
48" or 42" Waste Hose length
1/4" or 3/8" Inflow Hose
Standard Hose or G-Hose


An independent laboratory
study was conducted to det
bio-compatibility on the spe
and hosing. The laboratory
looking for harmful: chem


Designed for easier insertion and greater client comfort. Obturator is also designed with a longer thumb tab for ease in handling.

Lubricant packet (net wt. 5g) included.

An independent laboratory study was conducted to determine bio-compatibility on the speculums and hosing. The laboratory was looking for harmful: chemicals and toxicity that would be leached from the plastics itself. There was no findings of resins or toxicity found that could be derived by contact from this material. It was given a Superior Medical Rating.

BULB STYLE (PATENT NO. US D452004S) STRAIGHT STYLE







3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

NEVER RUN OUT OF STOCK WITH OUR **AUTOSHIP PROGRAM**

According to several management studies, the single most important characteristic of an effective office is the ability to manage time and worries.

Ask yourself these questions:

- Are you running out of stock?
- 8 Do you have a hard time meeting deadlines?
- 8 Do you hate keeping up with your inventory?
- 8 Do you call to overnight items, because you're running short?
- © Do you worry about canceling appointments, because you do not have what you need?
- © Do you hastily call to order stock, pressed for time, only to remember something more after you already hung up the phone?
- Are you busy with patients and do not have the time to call in the orders?
- 18 Were you ever "out of business" because your stock hasn't arrived or you forgot to place that order?
- 8 Do you check with everyone to make sure someone has placed your order?

If you answered YES to even one or more of these questions, you're sure to benefit from our popular, effective method to ship what you need, and have it there, when you need it.

The program I just introduced is our AUTO-SHIP Program.

- You choose what you want shipped.
- © When you want it shipped.
- How much you want shipped.



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

You now have peace of mind, NO stress, NEVER run out. Results: SAVE TIME and MONEY.

Use the enclosed "Auto Ship Program Agreement" in the pocket of this catalog to tailor your Auto-Shipment needs and drop it in the mail, e-mail it, or give us a call and we will help you get started.

Some suggestions along with your speculum order are, ordering a case of Underpads every other order, or adding Filters and gloves as needed. We can also take a look at your ordering history, which can help you in your decision.

We want to show you how to eliminate an unnecessary interruptions, so we suggest that you save time by making time to plan your Auto-Shipment right now, while you still have a mind to save time. Your business will be exploding!!

We Have It All Under Control.





3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Therapy Table & Accessories

Adjustable Electric Therapy Table





ADVANTAGES:

(Qualifies for Wheel Chair Accessible Tax Credit)

- Lowers to Wheel Chair Height
- Adjustment for Patient and Therapist Body Mechanics
- A Fully Equipped Bodywork Table
- Minimal Maintenance

FEATURES:

- Available in 30" wide x 6' long
- 10" Extension can be added to the table
- Two Designs (Flat or Tilt Back)
- Upholstery color is Buff
- Frame color is a textured ice white epoxy paint

SPECIFICATIONS:

- Top of the pad height at lowest point is 17"
- Top of the pad height at highest point is 30"
- Pad is 3" of med density foam
- Air operated foot controls
- Standard 115 VAC, 4 AMPS
- Standard 8 ft power cord



3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Light Weight Portable Table and Accessories





- Upholstery Standard Color Blue or Black (In Stock) 28 to 29 lbs
- Other Upholstery Colors By Special Order Only
- **Table Cart Carrier** Weight 9 lbs, but will hold 180 lbs
- 27" x 14" x 7" When Closed
- Enjoy the ease of rolling your table where ever you go.



- Pillows: Provide the correct level of support for prone, supine and side-lying positions
 - Wedge With Matching Vinyl Cover
 - ► **Round Bolster** 6" x 6" x 26"
 - ► Large Round Bolster 9" x 9" x 26"

 (Support under the knees, ankles or use for side posture support)









3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Therapist's Chair and Accessories

- Designed for the professional therapists
- Designed for comfort of sitting for long periods
- Seat in your choice of color
- Easy rolling on caster wheels



- Lightweight, yet sturdy
- "No-slip, no-tip" design
- Durable top, covered with black ribbed rubber and trimmed with grey molding



- Waterproof Table Cover
- Protects your table from moisture and grime
- No liquids seep into your table's foam



- Table Warmer with Thermostat Controls
- 2 styles Basic and Deluxe
- Keep clients warm, comfortable, and happy during treatment









3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Colon Therapy Charts






COLON MODEL



COLON THERAPY CHART





THE DIGESTIVE SYSTEM

BLUEPRINT FOR HEALTH




3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com






DISEASES OF DIGESTIVE SYSTEM

GASTROESOPHAGEAL DISORDERS






BOWEL REFLEX POINTS

(IBD) INFLAMMATORY BOWEL DIS



3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com



UNDERSTANDING IRRITABLE BOWEL SYNDROME (IBS)



STAINLESS STEEL VIEW TUBE BRUSH HOLDER


GOWNS


- **♦** Custom Made
- ♦ Small, Medium, Large, X-Large, and XX-Large
- **♦** Two different Lengths (To the Calves, or ankles)
- **♦** Fade-Resistant, Wrinkle-Free
- **♦** Cotton/Poly Seersucker Material
- **♦** Soft Elegant Pastel Stripes
- **♦** Velcro Closures





3145 SW 74th Terrace, Ocala, FL 34474
Order Lines 352-401-0303
www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com





THESE TWO GREAT NON-GLUTEN RECIPE BOOKS ARE A MUST FOR YOUR COOKBOOK LIBRARY! BY GWYNN PALMER





SPECULUM SAMPLE PACKET

SUPPLY PACKET



EXAM GLOVES



Attention Equipment and/or Speculum Buyers:

In accordance with Federal Regulations, it is a requirement that a letter be placed in the file of the equipment and/or purchasing speculums. When you purchase our equipment, we will need one for the equipment and one for the ordering of the disposable speculums. A sample letter has been enclosed for you to follow exactly with no variations for both in one script.

"Caution: Federal law restricts this device to sale by or on the order of a physician." 21CFR801.109

This regulation applies to all Class II Colon Therapy Equipment and Disposable Speculums by every Manufacturer both United States and Foreign Countries.

We also have a dual device listing on our disposables speculums for Class I devices, which has different package labeling. To ensure you are in compliance with the Federal Law, we will need to be advised of the class of device that you are using so that we may provide you with the proper package labeling.

Sample Letter

Letter Head with Doctor's Name, Address and Phone number

Clearwater Colon Hydrotherapy, Inc. 3145 S.W. 74^{th} Terrace Ocala, Florida 34474

Subject: Order to Purchase

Re: PPC-101 Colon Hydrotherapy Unit and/or Disposable Speculums

In Accordance with the Federal Code of Regulations 21 Food and Drugs, Parts 800-1299 as stated below.

"Caution: Federal law restricts this device to sale by or on the order of a physician." 21CFR801.109

This letter is to serve as a order for the following:

Person: _____Clinic: _____to purchase PPC-101 Colon Hydrotherapy Unit and/or Disposable Speculums/Kits.

Sincerely, Doctor's Signature Printed Name and License #



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

Product Information and Price List

ULTRA-VIOLET WATER TREATMENT SYSTEMS
Ultra-Violet Lamp (Replacement Bulbs 10000 Hour UV Bulb \$49.95
DDOEECCIONAL HWDDOTHED ADW MODELC
PROFESSIONAL HYDROTHERAPY MODELS
PPC-101 The Traveler
PPC-101 Clinical
PPC-101 Clinical with Specimen Collector
Options For Traveler:
Stainless Steel Cart with Sewage Tank
Stainless Steel Table (27 1/4" high)
DISPOSABLE SPECULUMS
Large & Small BULB Styles
Large, Midi & Small STRAIGHT Styles\$4.15
Kits with D/C Connectors
Kits Through the Auto-Shipment Program
Speculum Sample Pack - Standard\$19.75
Speculum Sample Pack - D/C \$20.75
(The Sample Packets Consists of: 5 Complete Kits - 1 Kit for Each Size and Style of
Speculums)
WATER FILTER REPLACEMENT CARTRIDGES
Granular Activated Carbon Cartridges (0.5 Micron Taste/Order/Chemical)
10 Inch
5 Inch
Fiber Sediment Cartridges (Rust/Sediment 1 Micron)
10 Inch
5 Inch
Carbon Filter (GAC Inline) For Other Manufacturer Units
10 Inch

Clearwater Colon Hydrotherapy, Inc. Product Catalog - Rev. 9.3 3-29-2023

Page 30



CLINICAL & PORTABLE COLON THERAPY TABLES AND ACCESSORIES

Light Weight Portable Table (28 to 29 Pounds)Call in for Updated PricingAdjustable Electric Therapy TableCall in for Updated PricingAdjustable Electric With Tilt OptionCall in for Updated PricingTable Cart Carrier\$99.00Plastic Sheet (Table Cover)\$28.95Table Cover Warmer (Basic)\$89.00Table Cover Warmer (Deluxe)\$120.95
Therapist's Chair\$425.00Pillow Wedge (with matching cover)\$180.00Bolster Full Round (with matching cover)\$65.95Bolster Jumbo Round (with matching cover)\$79.00Step Stool (Medical Grade Construction)\$69.00Toilet Step\$69.00
COLON THERAPY SUPPLIES Supply Packet (10% Discount Applied)
Order Neutralizer/Cleaner, Pack of 20 Sanitation End Connectors, Stainless Steel Cleaner, Surface Disinfectant, 2 Sediment Filters, and 2 Carbon Filters) Surface Disinfectant Decontaminate Cleaner
Sanitation End Connectors (Pack of 20)
View Tube Brushes (Pack of two 24" Length) \$24.00 Stainless Steel Brush Holder \$59.95
View Tube for PPC-101 Clinical Model (Quantity 1)

Clearwater Colon Hydrotherapy, Inc. Product Catalog - Rev. 9.3 3-29-2023

Page 31



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303

www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

•	Bulb for View Tube Assembly (Quantity 1) 9000 hr Natural Lighting \$18.00 Gloves - Powder Free, Nitrile Gloves
Unde	rpads (Large Case of 150)
Dispo	osable Footies
	Large - 42"
	Large - (Custom 51" Size to Ankle)
	cant (Tube)
	mmended Reading for Certification. These are Excellent References to have in library. Items Marked with an * are Suggested.
	BOOKS Office Safety Manual (Set to OSHA Standards for Colon Therapy
your :	BOOKS Office Safety Manual (Set to OSHA Standards for Colon Therapy Author: Mary Ruth Baker
your	BOOKS Office Safety Manual (Set to OSHA Standards for Colon Therapy



3145 SW 74th Terrace, Ocala, FL 34474 Order Lines 352-401-0303 www.colonhydrotherapy.com ~ E-Mail: sales@colonhydrotherapy.com

CD's

CD S	
Medicine Woman (CD) Part I	\$15.95
Medicine Woman (CD) Part II	
Medicine Woman (CD) Part III	
Wheat Belly - (Set of 6 CD's)	
Why People Don't Heal and How They Can - (Set of 6 CD's)	
MODELS	
Colon Model	\$82.80
CHARTS	
Colon Therapy (Chart of the Colon)	\$23.99
Digestive System	\$23.99
Blueprint for Health Chart	
Diseases of the Digestive System	\$23.99
Gastroesophageal Disorders	
Inflammatory Bowel Disease (IBD)	
Understanding Irritable Bowel Syndrome (IBS)	
Quick Check List - Morning Set-up, Therapy & Disinfecting the Unit	After Use
For Clinical	
For Traveler	
1 Of Traveler	

PLEASE NOTE: Prices are subject to change without notice and do not include shipping charges.